ICT

Operating Instructions

ICE COMPRESSION DEVICE / KNEE - LEFT OR RIGHT KNEE (10010 001-00 / 10110 040-00)

2-part system provides optimal cold therapy for knee pain and inflammation from strains/sprains, ACL / MCL / Meniscus strains or tears, patella tendonitis (jumper's knee), iliotibial band friction (runner's knee). Customizable fit for maximum compression and comfort*.

ICE COMPRESSION DEVICE / BACK* - LOWER AND MIDDLE BACK (10040 001-00)

2-part system that provides optimal cold therapy for back pain associated with any of the following: soreness, stiffness, spasms, overuse, strain, arthritis, SI joint pain, sciatica. Customizable fit for maximum compression and comfort.

ICE COMPRESSION DEVICE / UTILITY - ANKLE, ELBOW, WRIST, SHIN (10030 001-00 / 10130 040-00)

2-part system provides optimal cold therapy for high or low ankle sprains, elbow tendonitis (tennis elbow), wrist sprains, shin splints, plantar fasciitis, bruised feet or hands, achilles strain / tendonitis, calf strains, bruised forearms and biceps. Customizable fit for maximum compression and comfort*.

ICE COMPRESSION DEVICE / SHOULDER - RIGHT OR LEFT SHOULDER / SOLD SEPARATE (10021 001-00 / 10022 001-00 / 10121 040-00 / 10122 040-00)

2-part system that provides optimal cold therapy for rotator cuff tendonitis / tears, AC joint pain, strained tendons, bursitis, shoulder impingements and overuse of the shoulder. Customizable fit for maximum compression and comfort*.

ICE COMPRESSION DEVICE / EXTENDED SHOULDER* - RIGHT OR LEFT SHOULDER / SOLD SEPARATE (10051 001-00 / 10052 001-00)

2-part system that provides optimal cold therapy for rotator cuff tendonitis / tears, AC joint pain, strained tendons, bursitis, shoulder impingements and overuse of the shoulder. Customizable fit for maximum compression and comfort. Includes Oversized Ice Cell for double area coverage.

*Hyperice Pro only

USAGE INSTRUCTIONS

- 1) Fill the Ice Cell 2/3 with ice or HYPERICE FUEL™ – do not add water. Crushed ice or small ice nuggets work best with all HYPERICE® products, however ice cubes may be substituted as needed.
- 2) Place the Ice Cell on a flat surface and flatten out - removing as much air from the Ice Cell as possible before securing the cap.
- 3) Insert the Ice Cell into the Compression Wrap.
- 4) Place on the body and apply the compression straps to the desired level of compression.
- 5) Lightly press the Air Release Valve to remove air from the Ice Cell. For added compression, re-tighten the compression straps after the air has been released.
- 6) Leave the device on for 15 – 20 minutes. Press the Air Release Valve once every 5 minutes to release any air that accumulates during the treatment to sustain optimal compression.

CARE INSTRUCTIONS

- 1) After use, rinse the Ice Cell inside and out with fresh water and air dry.
- 2) After the Ice Cell is dry, store in a cool, dry place.
- 3) Do not expose the Ice Cell to direct sunlight and/or heat, as the UV rays can damage the medical-grade material.
- 4) Do not fill the Ice Cell with anything other than ice.
- 5) Do not machine wash or iron.

CAUTION

- 1) This product contains natural rubber latex which may cause allergic reactions.
- 2) Do not leave any HYPERICE® device on the body for more than 25 minutes at a time.

For more detailed directions on how to use your HYPERICE® device, please visit the PRODUCTS section of our website at www.hyperice.com. There you will find an instructional video that clearly illustrates how to use your HYPERICE® device.

1-YEAR WARRANTY

If within one year from the date of purchase, this product fails due to a defect in material or workmanship, Hyperice Inc., Hyperice SAS and/or the international distributor, dealer or retailer from whom the purchase was made, will repair or replace the product, or necessary components, free of charge.

THIS WARRANTY EXCLUDES:

- (a) damage caused by accident, abuse, mishandling, or transport;
- (b) units subjected to unauthorized repair;
- (c) units not used in accordance with Hyperice care instructions;
- (d) damage exceeding the cost of the product;
- (e) deterioration of the delivered product resulting from abnormal storage and/ or safeguarding conditions on the client's premises;
- (f) failure to provide the dated proof of purchase.

Some states and countries do not allow a limitation of damages, so the foregoing limitation may not apply to you. This warranty guarantees specific legal rights, and but other rights may vary from country to country and from state to state. This warranty service is available through the operating policies and procedures of Hyperice Inc., Hyperice SAS, and all Hyperice international distributors. If you are a client, please contact the appropriate Hyperice business entity, dealer or retailer, from whom you directly purchased the product for warranty and return authorization procedures. This warranty is initiated and executed by Hyperice Inc. and Hyperice SAS.

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Copyright 2021 © Hyperice, Inc. All rights reserved.
Hyperice, Inc.